


Who Dunnit?


Purpose: To compare the fibers taken from suspects to a sample fiber found at the crime scene to determine WHO DUNNIT!

Materials:



Procedure: Obtain a "suspect" fiber sample and prepare a WET MOUNT.

- Why do we prepare a wet mount? _____



- Repeat steps for all suspects.
- Compare microscope drawings to find out WHO DUNNIT.

Observations: ****Note: Observe under HIGH power****

Crime Scene	Power _____ x
	Description: _____

Suspect _____	Power _____ x
	Description: _____

Suspect _____	Power _____ x
	Description: _____

Suspect _____ Power _____ x

Description: _____

Suspect _____ Power _____ x

Description: _____

Conclusions & Questions:

1) Could you figure out WHO DUNNIT without the microscope? Explain. (2 marks) _____

2) Based on your observations and drawings, WHO DUNNIT? Give at least 2 reasons how do you know. (3 marks)

